

### **IN THE SPECIFICATION**

Amend the paragraph beginning on page 8 at line 3, as follows:

Two holes are drilled at the end of the bone at a spacing corresponding to the width of the implant as measured by the spacing of the legs 4 of the implant device thereof. The legs 4 of the implant device are impacted longitudinally into the drilled holes entering and aligning to the medullary canal. The distal ends of the legs 4 are tapered as shown in the drawings to facilitate the impacting of the legs in the bone. The fracture site is closed and the implant device is firmly seated and secured with the bone screw and washer to the bone at one end of the implant device. Compression at the fracture is achieved by turning the cam between the washer and the U-shaped bend of the implant device to effect further compression whereafter the screw is fully tightened and the washer is seated and then the cam is removed. In lieu of the cam, the tension force in the wire element can be produced by the surgeon applying pressure to the U-shaped bend portion 11 and then tightening the bone screw 16 while the wire is under tension.